

Im Hofstätt 11
D-72359 Dotternhausen
Tel +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
e-mail info@proactiv-gmbh.de
www.proactiv-gmbh.de



TRAVELER

Active easy running rigid frame wheelchair

Maximum user weight limit 120 kg

A9030100000 € 4,290.00

Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____

Order date: _____
according to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: **same as company address**
 Other (see below)

Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:

Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.O.B: _____
Paying organisation: _____

Invoicing address: **same as company address**
 other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Basic size

Body size:


Body height: _____ cm

Weight: _____ cm

Lower leg length: _____ cm


(anatomic dimension: hollow of the knee to the foot sole)

Indication: _____



This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Wheelchair dimensions:



Seat width SB	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	Optional with long frame:			
	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48	<input type="checkbox"/> 49	<input type="checkbox"/> 50	<input type="checkbox"/> 51	<input type="checkbox"/> 52
Seat depth ST	Short frame					Long frame						
	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44		
Back height RH	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29		
	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39		
	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48			
Lower leg length USL	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44				
	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48	<input type="checkbox"/> 49	<input type="checkbox"/> 50	<input type="checkbox"/> 51					
Seat height back SHh (for 24" drive wheels)	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input checked="" type="checkbox"/> 43					
Seat height front SHv (for 5" caster wheel)	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input checked="" type="checkbox"/> 47	<input type="checkbox"/> 48	<input type="checkbox"/> 49	<input type="checkbox"/> 50	<input type="checkbox"/> 51	<input type="checkbox"/> 52			
Tipping point KP (front side rear tube to axis middle drive wheel)	<input type="checkbox"/> 4	<input type="checkbox"/> 5	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input checked="" type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10	<input type="checkbox"/> 11	<input type="checkbox"/> 12	<input type="checkbox"/> 13		

Please select the required dimensions (all dimensions are in cm).

The underlined values are the standard settings.

The values SHh and SHv may deviate from the specified values depending on the wheel size and wheel camber!

Special dimensions possible upon request (except seat width).

Frame colour - one-colour design

Information on current colour designs and example images at: www.proactiv-gmbh.de/farben.php

Frame colour: Selection from colour palette

Frame colour will be specified with order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

Powder colours

White RAL 9016

gloss A9000100170-023

matt A9000100170-123

€ 42.10

Brilliant silver

gloss A9000100170-011

Anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101

€ 42.10

Black RAL 9005

gloss A9000100170-003

matt A9000100170-019

€ 42.10

Pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122

€ 42.10

Pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137

€ 42.10

Pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136

€ 42.10

Lemon green

gloss A9000100170-005

matt A9000100170-105

€ 42.10

Yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146

€ 42.10

Red RAL 3001

gloss A9000100170-047

matt A9000100170-147

€ 42.10

Lagoon blue

gloss A9000100170-052

matt A9000100170-152

€ 42.10

Athletic grey

gloss A9000100170-050

matt A9000100170-150

€ 42.10

Muddy grey

gloss A9000100170-051

matt A9000100170-151

€ 42.10

Powder colours „Dormant“

Deep black effect silver

gloss A9000100170-007

Sky blue

gloss A9000100170-017

Navy blue

gloss A9000100170-045

matt A9000100170-145

€ 42.10

Electric green

gloss A9000100170-015

matt A9000100170-115

€ 42.10

Violet

gloss A9000100170-044

Fire red

gloss A9000100170-018

matt A9000100170-118

€ 42.10

Frame colour - One-colour design special colours “liquid coating”

€ 339.40

Storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

Coffee

gloss A9000100170-520

matt A9000100170-521

Prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

Pearl (mother of pearl-white)

gloss A9000100170-524

matt A9000100170-525

Berry

gloss A9000100170-526

matt A9000100170-527

Blue ocean

gloss A9000100170-522

matt A9000100170-523

Emerald

gloss A9000100170-514

matt A9000100170-515

Mint

gloss A9000100170-528

matt A9000100170-529

Chili

gloss A9000100170-516

matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame) € 260.70

Please select a color from below:

Base colour (on the rear frame)

White RAL 9016

A9000100170-200

Brilliant silver

A9000100170-201

Anthracite RAL 7016

A9000100170-202

Black RAL 9005

A9000100170-203

Pastel blue RAL 5024

A9000100170-206

Lemon green

A9000100170-207

Yellow RAL 1023

A9000100170-208

Lagoon blue

A9000100170-212

Athletic grey

A9000100170-210

Muddy grey

A9000100170-211

Design colour (on the front frame)

Black RAL 9005

gloss A9000100170-259

matt A9000100170-260

Lemon green

gloss A9000100170-250

matt A9000100170-251

Chili

gloss A9000100170-253

matt A9000100170-254

Pure orange RAL 2004

gloss A9000100170-255

matt A9000100170-256

Coffee

gloss A9000100170-278

matt A9000100170-279

Mint

gloss A9000100170-286

matt A9000100170-287

Emerald

gloss A9000100170-274

matt A9000100170-275

Storm (ice blue)

gloss A9000100170-264

matt A9000100170-265

Prairie (rose gold)

gloss A9000100170-266

matt A9000100170-267

Pearl (mother of pearl-white)

gloss A9000100170-282

matt A9000100170-283

Berry

gloss A9000100170-284

matt A9000100170-285

Blue ocean

gloss A9000100170-280

matt A9000100170-281

Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

Drive wheel: aluminum rim, hub and spokes

Black A9000100171-003

Caster fork and bearing block

Silver A9000100173-011

Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: Aluminium rim and design inlet¹

Frame colour A9000100171-999 € 242.30

Select from one-colour design (different from the frame color) A9000100171-998 € 280.50

Colour: _____

A9000100170-

Caster wheel: 5" aluminium rim² (3-spoke design)

Frame colour A9000100174-999 € 73.00

Select from one-colour design (different from the frame color) A9000100174-998 € 84.60

Colour: _____

A9000100170-

Special colour requests (frame or attachment parts)

Frame colour for single purchasing Selection from RAL colour palette

A9000100170-100 € 402.30

Gloss level

gloss

matt

colour name

RAL number

Note: attachment parts are **not** possible in design colors

¹ Only in conjunction with Liteline hub (with drive wheel aluminum rim with 18 or 36 spokes, single rim)

² Only in conjunction with caster wheel 5" aluminum rim 3-spoke design good shock absorption A900050121*, as well as with suspension fork 5-spoke design A9000500160*

Frames

Straight front end


A9030100010-K*



Front end V-shape

A9030100020-K* € 294.80

Note: Footrest area with V-shape
x = SB - (0.25 x USL + 6 cm)



Custom-made seat width dimension 49 - 52 cm

A9030100060 € 297.00

Frame angle

Knee angle 105°

A9000100050-001

Knee angle 95°

A9000100050-002 € 329.10

Rear frame angle 105°

A9030100100-001
KP from 6.5 cm to 13 cm

Rear frame angle 95°

A9030100100-002 € 125.40
KP from 4 cm to 10 cm

Push handles

Back tube with integrated handles

A9000401051-K € 126.60
straight



Aluminum push handles fixed in back tube

A9000401050 € 126.60
from RH 23 cm



Horizontal push handles screwed into the back tube

A9000401076 € 234.80
can be unscrewed



Safety push handles height adjustable

A9000401060* € 210.10

Mount in the back tube
from RH 30 cm with adaptable back
from RH 35 cm for back tube with lordosis
support
from RH 43 cm with back shell



Safety push handles back-positioned length 45 cm


A9030401071 € 396.40

continuous height adjustment,
Safety push handle holder mounted on the
back cross bar
Up to SB 42 cm or body weight 95 kg
Not in combination with swing-away armrest



Equipment without push handles

A9000401000



Backrest

Folding backrest

A9000400011
7-way angle adjustable,
with integrated plastic-coated aluminum side
parts



Back tube with lordosis support

A9000400005-K € 111.10

- 1 cm
- 2 cm
- 3 cm



Back equipment

If not specified, the "Adaptable back with back padding" will be delivered.

Adaptable back with back padding

A9000802010-K
Cover: Black, breathable, open surface
Cushion thickness: 16 mm
Equipment with back cross bar



Adaptable back with back padding „light“

A9000802320-K001 € 44.90
Back padding with extra thin padding for easy, space-saving folding of the wheelchair.
Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment with back cross bar

Back rise + 3 - 5 cm

A9000802100 € 151.60
for the nominal back height dimension
Only in combination with height-adjustable push handles



Equipment without back upholstery and back padding or back system

A9000802011 € - 53.00
Equipment with back cross bar

For other back systems for example ergonomic back shell, see separate order sheet: "PRO ACTIV seat and back systems"



You will find our PRO ACTIV transport restraint systems on our separate "Transport restraint systems" order sheet




Seat upholstery

If not specified, the "Open belt system" will be delivered.

Open belt system

A9000801011-K



Closed belt system

A9000801020-K € 89.70





Equipment without seat upholstery

A9000801050



Seat cushion



Seat cushion in rectangular contour Inlet: Standard

A9000805180-K* € 111.90

Seat cushion with radius contour Inlet: Standard

A9000800300-K* € 155.40

Seat cushion height:

3 cm 5 cm
 4 cm 6 cm

Seat cushion in rectangular contour Inlet: Anti-Decubitus

up to 30x30 cm A9000805204-K* € 290.00

>30x30 up to 50x50 cm A9000805200-K* € 422.00

Seat cushion with radius contour Inlet: Anti-Decubitus

up to 30x30 cm A9000800311-K* € 342.80

>30x30 up to 50x50 cm A9000800310-K* € 465.90

Seat cushion height:

4 cm
 6 cm

Cover: Black, breathable, open surface Black, breathable, closed surface

If no specification is made, the cover material "Black, breathable, open surface" will be delivered.

 Seat cushion is existing**You can also find the seat cushions on our separate order sheet: "PRO ACTIV seat and backrest systems"****Clothing guard****Aluminum clothing guard removable via lock function**

Adjustable following wheel arch

Only in combination with folding backrest



in anthracite matt

 A9000403210-K € 236.10

in frame colour (basic colour with two-colour design)

 A9000403211-K € 302.00

in design colour (only possible with frame colour in two-colour design)

 A9000403212-K € 335.00**Carbon clothing guard removable via lock function** A9000400012-K € 547.30

Adjustable following wheel arch

Only in combination with folding backrest

**Armrests****Armrest, height adjustable, integrated in the clothing guard**

Including aluminum clothing guard

with arm padding 25 cm long

Removable via lock function

Only in combination with folding backrest



in anthracite matt

 A9000400008-K € 553.00

in frame colour (basic colour with two-colour design)

 A9000400009-K € 603.80

in design colour (only possible with frame colour in two-colour design)

 A9000400010-K € 629.20**Armrest, height adjustable, integrated in the clothing guard** A9000400018-K € 836.90

Including carbon clothing guard

with arm padding 25 cm long

Removable via lock function



Only in combination with folding backrest

Axle plate

Standard equipment with 12 mm sockets/axles. Further axle plates on request.

2-slot axle plate A9000201000*Seat height back = approx. 37 - 43 cm
(for 24" wheel)**2-slot extended axle plate** A9000201010* € 169.60

Wheelbase extension by 7.5 cm

**Extended bike plate** A9000201030* € 306.00+ 4° camber integrated,
Bike track can be set separately

Wheel camber

*Total width with tight-mounted handrim.

1° drive wheel camber

A9000202010*
Total width = approx. SB + 20 cm*

4° drive wheel camber

A9000202030*
Total width = approx. SB + 26 cm*

6° drive wheel camber

A9000202040*
Total width = approx. SB + 30 cm*

E-Adaptations

Preparation for the "e-motion" / „e-motion DuoDrive“ mounting

A9000201200 € 488.30

Equipment: e-fix torque support, 1/2" axle and socket length 50.8 mm (2")
PRO ACTIV and Alber recommend anti-tipping supports are installed in pairs. Please choose "Adaptations" under the selection.

Preparation for the "e-fix" mounting

A9000201210 € 488.30

Equipment: e-fix torque support, 1/2" axle and socket length 50.8 mm (2")
PRO ACTIV and Alber recommend anti-tipping supports are installed in pairs. Please choose "Adaptations" under the selection.

Brake system

If not specified, the "Brake system push to lock with standard brake handle" will be delivered.

Recommendation: For KP 7 cm or lower, we recommend the „Brake system push to lock with brake handle 20° outwards“.

Brake system push to lock Standard brake handle

A9000701010-001



Brake system push to lock Long brake handle

A9000701010-002



Brake system push to lock Low brake handle

A9000701010-003



Brake system push to lock Brake handle 20° outwards

Standard brake handle
A9000700002

Long brake handle
A9000701010-004



Brake system push to lock Brake handle to flip-down


A9000706008 € 154.00
Length: 9 cm



Drum brake

A9000702011-K* € 460.20

Operating from the push handle,
Only with 2-slot axle plates
(Total width (see "Camber" selection) + 3 cm
wheel plate cannot be mounted on the inside).
Please note:
- Only in combination with drive wheel
aluminum 36-spoke rim wheel
- Not in combination with one-hand drive,
e-fix or e-motion



Footrest

If not specified, the "continuous footrest" will be delivered.

Continuous footrest

A9030601100-K
with tubular resting bracket, foldable,
removable
When folding the wheelchair, the footrest also
folds completely from USL 37 cm.
With foot support area x < 22 cm the foot
bracket is designed to be parallel.
If the footrest is wanted with a "rubber coated
footplate", please select this under "Footrest
accessories"



Footrest folding up to one side

A9030601200-K* € 240.10

with continuous tubular resting bracket,
height adjustable +/- 3 cm,
can be pivoted outwards, removable
V-shape from SB 40 cm
If the footrest is wanted with a "rubber coated
footplate", please select this under "Footrest
accessories"



Footrest folds up to the rear with spring locking mechanism

A9030602005-K* € 279.00

with continuous tubular resting bracket
for KW 105° short frame from USL 37 cm, long
frame from USL 39 cm
for KW 95° short frame from USL 36 cm, long
frame from USL 38 cm
If the footrest is wanted with a "rubber coated
footplate", please select this under "Footrest
accessories"






<p>Divided footrest</p> <p><input type="checkbox"/> A9030603100-K € 392.90</p> <p>with rubber coated footplate, angle-adjustable, can be pivoted outwards and folded up to the side, removable V-shape from SB 40 cm</p> 	<p>"Swing Away" footrest</p> <p><input type="checkbox"/> A9016604100-K € 392.90</p> <p>with rubber coated footplate, separated in the middle, pivoted to the side, can be folded up, height adjustable, removable, angle-adjustable</p> 	
---	---	--


Footrest accessories

<p>Rubber coated footplate</p> <p><input type="checkbox"/> A9000610000-K € 122.20</p> <p>bolted, angle-adjustable Depth: 14 cm</p> 	<p>Footrest side stop</p> <p><input type="checkbox"/> A9000610100 € 58.00</p> <p>adjustable only in combination with rubber coated footplate</p> 	
---	--	--

Adaptations

<p>Anti-tipping support</p> <p><input type="checkbox"/> right side A9000300110-K* € 164.20</p> <p><input type="checkbox"/> left side A9000300111-K* € 164.20</p> <p>adjustable, removable, can be pivoted, with tension spring and universal adapter</p> 	<p>Transit wheels</p> <p><input type="checkbox"/> A9000301110-K* € 320.00</p> <p>with universal adapter, removable not in combination with anti-tipping support</p> 	<p>Tilting clamp</p> <p><input type="checkbox"/> right side A9000302110* € 136.00</p> <p><input type="checkbox"/> left side A9000302111* € 136.00</p> <p>with universal adapter and Quick pin Not in combination with anti-tipping support, right side or left side Colour of the mount: <input type="checkbox"/> black <input type="checkbox"/> silver If not specified, the mount will be delivered in black.</p> 
--	---	---

3" caster wheel

<p>3" safety light-up wheel (80 mm)</p> <p><input type="checkbox"/> A9000501110 € 68.70</p> <p>puncture-proof</p> 	
--	--

4" caster wheel

<p>4" caster wheel „Standard“ solid rubber black</p> <p><input checked="" type="checkbox"/> A9000501011-004 (82A), puncture-proof</p> 	<p>4" Skate wheel (100 mm)</p> <p><input type="checkbox"/> A9000501120 € 68.70</p> <p>puncture-proof</p> 	<p>4" safety light-up wheel (100 mm)</p> <p><input type="checkbox"/> A9000501220 € 92.00</p> <p>puncture-proof</p> 
--	---	--

4" Frog Legs caster wheel (100 mm)

A9000501230 € 206.30

with replaceable rubber 6 mm, puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



5" caster wheel

5" caster wheel „Standard“ PU black

A9000501036-005
(82A), puncture-proof



5" caster wheel „Standard“ PU grey

A9000501031-005
(82A), puncture-proof



5" safety light-up wheel (125 mm)

A9000501221 € 92.00
puncture-proof



5" aluminum rim 3-spoke design “good shock absorption”

Silver
A9000501210 € 174.60

Black
A9000501211 € 174.60

with replaceable rubber 6 mm medium hard,
puncture-proof
Incl. hair guard spacer (3 mm)



5" Frog Legs caster wheel (125 mm)

A9000501232 € 206.30

with replaceable rubber 6 mm, puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



6" caster wheel

6" caster wheel „Standard“ PU grey

A9000501010-006
(82A), puncture-proof



6" caster wheel “pneumatic tyre”

A9000501050 € 55.50
SH front min. 42 cm



6" caster wheel “Flexel System”

A9000501070 € 55.50
puncture-proof



Caster fork

If not specified, the “standard caster fork with screw axle” is supplied.

Caster fork Standard

A9000500100-K001
with screwed axle

A9000500100-K002 € 162.50
with quick-release axle



Caster fork Heavy Duty for high loads

Extra weight approx. 100 - 150 gr.

A9000510111-K001 € 52.40
with screwed axle

A9000510111-K002 € 212.50
with quick-release axle



Suspension fork

inkl. 5" aluminum rim “good shock absorption”

A9000500160-K002 € 598.90
with quick-release axle



Caster fork spare wheel „Standard“

A9000500300-K € 288.10

Caster fork Standard with quick-release axles



Caster fork spare wheel „Heavy Duty“

A9000500301-K € 328.10

Caster fork Heavy Duty with quick-release axles




The spare caster fork does not have a caster wheel.

In case a caster wheel is requested, please put a cross on a second selection in the category "Caster wheels".

Track distance

Caster wheel track extension



A9000500400* € 57.40

10 mm per side

A9000500401* € 57.40

15 mm per side

A9000500410* € 57.40

20 mm per side

Drive wheels

Information about other rim types, designs and availability can be found at: www.proactiv-gmbh.de.
If not specified, the "drive wheel 36 spokes aluminum rim" will be delivered.

Drive wheel 18-spoke aluminum rim

24" (ETRTO 540)
A9000203122

with Liteline hub, radially spoked
Not in combination with:

- Puncture-proof tyres
 - One-hand drive
 - Drum brake
 - Transport restraint connection
 - Handbike
- up to 70 kg body weight:

Gewicht	Sturz				
	1°	4°	6°	8°	10°
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



Drive wheel 36-spoke aluminum rim

24" (ETRTO 540)
A9000203000-010

with Liteline hub (left image),
radially spoked
up to 120 kg body weight

24" (ETRTO 540)
A9000203000-020

with aluminum front hub (right image),
radially spoked
over 120 kg body weight



Drive wheel double rim 36-spokes

24" (ETRTO 540)
A9000203100 € 153.20

25" (ETRTO 559)
A9000203200 € 153.20

26" (ETRTO 590)
A9000203300 € 153.20

with aluminum front hub,
radially spoked
Use with higher loads and body weight
> 120 kg



Spinergy wheel "Spox Sport" with 24 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)
A9000203550 € 761.30

25" (ETRTO 559)
A9000201559 € 761.30

26" (ETRTO 590)
A9000203555 € 761.30



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim and hub in black, spokes in white

24" (ETRTO 540)
A9000203554 € 929.30

Recommendation: up to 70 kg body weight,
not in combination with handbike/traction
device use.



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim, hub and spokes in black


24" (ETRTO 540)
A9000203552 € 929.30

Recommendation: up to 70 kg body weight,
not in combination with handbike/traction
device use.



Handrims

x = dimension handrim attachment



Handrim attachment
general:

- x = 11 mm
- x = 19 mm

Handrim attachment
Curve + Gekko handrim:

- x = 10 mm
- x = 18 mm

If not specified, the handrim will be mounted wide (19 mm / 18 mm).

Stainless steel Ø 20 mm

- 24" handrim for 24" wheel and 25" wheel A9000205010-004
- 26" handrim for 26" wheel A9000205010-006



Aluminum Ø 19 mm silver anodised

- 24" handrim for 24" wheel and 25" wheel A9000205019-004



Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 24" handrim for 24" wheel and 25" wheel **€ 66.30** A9000205020-004
- 26" handrim for 26" wheel **€ 66.30** A9000205020-006



Aluminum Ø 19 mm polished and hard-anodised

- 24" handrim for 24" wheel and 25" wheel **€ 160.10** A9000205030-004
- 26" handrim for 26" wheel **€ 160.10** A9000205030-006




Titanium Ø 19 mm

- 24" handrim for 24" wheel and 25" wheel **€ 325.80** A9000205040-004
- 26" handrim for 26" wheel **€ 325.80** A9000205040-006



Silicone handrim cover

- for 24" handrim **€ 226.00** A9000205100-004
 - for 26" handrim **€ 226.00** A9000205100-006
- Surface without structure, compatible with handrim Ø 19 mm or Ø 20 mm.



Ultra-Grip „very grippy“

- 24" handrim for 24" wheel and 25" wheel **€ 390.70** A9000205070-004
- 26" handrim for 26" wheel **€ 390.70** A9000205070-006



Ergo-Grip

- 24" handrim for 24" wheel and 25" wheel **€ 396.80** A9000205080-004
 - 26" handrim for 26" wheel **€ 396.80** A9000205080-006
- Attachment fitting to the rim



Curve black anodised, ergonomic handrim profile

- 24" handrim for 24" wheel and 25" wheel **€ 389.10** H 24 mm A9000205130-004
- 24" handrim for 24" wheel and 25" wheel **€ 389.10** H 31 mm, size L A9000205140-004
- 26" handrim for 26" wheel **€ 389.10** H 31 mm, size L A9000205140-006



<p>Curve silver anodised ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel H 31 mm, size L € 394.30 A9000205160-004</p> 	<p>Curve powder-coated, ergonomic handrim profile</p> <p><input type="checkbox"/> = Frame colour (gloss surface) <input type="checkbox"/> = Black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel H 24 mm € 454.00 A9000205135-004</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel H 31 mm, size L € 454.00 A9000205145-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel H 31 mm, size L € 454.00 A9000205145-006</p> 	<p>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel H 31 mm, size L € 505.50 A9000205150-004</p> 
<p>Quadro with Tetra-Grip coating</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 517.50 A9000205170-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel € 517.50 A9000205170-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 517.50 A9000205170-006</p> 	<p>Gekko black anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 501.20 A9000205175-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel € 501.20 A9000205175-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 501.20 A9000205175-006</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile</p> 	<p>Gekko silver anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 513.20 A9000205185-004</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile</p> 
<p>Gekko S black anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 501.20 A9000205195-004</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm</p> 	<p>Equipment without handrims</p> <p><input checked="" type="checkbox"/> A9000205990</p>	

Spoke guard

Finger- and spoke guard "transparent"

for 24" wheel
A9000206300-004 € 201.10

for 26" wheel
A9000206300-006 € 201.10

Not in combination with Spinergy drive wheels.



Tyres

If not specified, the "Schwalbe RightRun in black" tyres will be delivered.

Schwalbe „RightRun“ High-pressure tyres

24x1" (540-25)
A9000204001-004

Colour: Black
7.5 bar



Schwalbe „RightRun“ High-pressure tyres

24x1" (540-25)
A9000204000-004

25x1" (559-25)
A9000204000-005

26x1" (590-25)
A9000204000-006

Colour: Grey, black wall
7.5 bar



Schwalbe „DownTown“ Profiled tyres

24x1" (540-25)
A9000204010-004

7.5 bar



Schwalbe „Marathon Plus Evolution“ with "smart guard" puncture protection

24x1" (540-25)
A9000204030-014 € 47.80

25x1" (559-25)
A9000204030-015 € 47.80

26x1" (590-25)
A9000204030-016 € 47.80

Colour: Black
10 bar



Schwalbe „Marathon Plus“ with "smart guard" puncture protection

24x1" (540-25)
A9000204030-004 € 47.80

Colour: Grey, black wall
10 bar



Schwalbe „One“

24x1" (540-25)
A9000204123-004 € 90.30

Colour: Black
9 bar



Puncture-proof "solid rubber tyres"

24x1" (540-25)
A9000204201-004 € 119.10

Colour: Black
Not in combination with the 18 spoke drive wheel, Spinergy and one-hand drive



Accessories

24" MTB spare wheel set

with double rim (ETRTO 507), aluminum front hub, quick-release axle and stainless steel handrim
Mountain bike tyres, selectively:

A9000203410 € 682.90

24x1.75" Schwalbe "Marathon Plus"
Brake setting compatible with 24x1" wheel



A9000203413 € 682.90

24x2.1" Schwalbe "SMART SAM Sport"
Brake setting not compatible with 24x1" wheel



NOTE: We charge a surcharge for the use of drum brakes

A9000200004 € 221.10

Through axle with assistance, Large push button

A9000200096 € 102.70

Interlock dimension: 122 / 110 mm
Axle diameter: 12 mm
Not in conjunction with lockable drive wheel axle clip "Tetra Clip",
not in combination with one-hand drive



Lockable drive wheel axle clip "Tetra Clip"

A9000206400 € 161.60

Not in combination with Spinerger drive wheels,
one-hand drive, drum brake or Liteline hub



Hand air compressor with pressure gauge

A9000901100 € 103.20

up to 14 bar



Length-adjustable safety belt with buckle

A9000902100 € 78.30




Adjustable safety belt with hook and fastener

A9000902110 € 62.60



Crutch holder with flexible nylon cup

A9000902310* € 133.40



SIGMA lighting set

A9000903270 € 104.80

Consisting of front light and rear lights,
chargeable via Micro USB socket,
incl. Micro USB charging cable
Lighting set compliant with road traffic regulations



Transport bag for Traveler-frame

A9030903200 € 237.00

Only in conjunction with caster forks with quick-release axles



Transport bag for drive wheels

for 24" wheel
A9000903530-004 € 75.50
Ø 76 cm x 2 cm

for 26" wheel
A9000903530-006 € 75.50
Ø 83 cm x 2 cm



Locking strap for Traveler

A9030902200 € 19.50

Further accessories (such as handbikes, steering attachments, etc.) see separate order form or www.proactiv-gmbh.de.

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 16 - 2024-06-01 / The prices are plus statutory value-added tax and shipping costs = Included in the basic price = Position selected
Order sheet TRAVELER Rev. 07 UK