



STRONG IN USE, FLEXIBLE IN HEIGHT.

The **LIFT solid** is a real "workhorse". Thanks to the infinitely variable seat height adjustment, either mechanically or electrically, it opens up flexible activities in a wide variety of professions and the use of conventional furniture and kitchen equipment. The LIFT solid shows its true strength in demanding environments: robust, resistant and ready for greater challenges.



#Independence

With the LIFT solid, you can master your everyday life and your professional life in a self-determined way. It gives you the freedom to act independently in almost any situation.



#AtEyeLevel

The infinitely variable seat height adjustment raises the seat by up to 30 cm. In this position, the wheelchair user is at eye level with a person standing about 1.75 m tall.



#Flexibility

Your LIFT solid adapts to you: The tipping point can be adjusted afterwards to meet your needs.



#Independence

Thanks to the foldable backrest and the tool-free removal of assemblies, you can easily transport your LIFT solid.





MADE FOR EVERY CHALLENGE.

Your **LIFT solid** will be manufactured individually for you, to suit your needs. You have a large selection of equipment – from footrests and brakes to push handles, back and seat systems as well as drive and caster wheels. This makes your LIFT activ exactly as you need it.

UNIQUE, SELECTABLE OPTIONS FROM PRO ACTIV:



The 7-way angle-adjustable and foldable backrest offers maximum flexibility for the ideal sitting position.



The removable and swivelling frame stem with split footrest offers plenty of space for your transfer.



The footrest, which can be folded towards the back, is easy to handle and flexible in everyday life.



With the Heavy Duty equipment, your LIFT solid is designed for high loads and is suitable for user weights up to 140 kg.

TECHNICAL SPECIFICATIONS

Seat width:	35 - 50 cm (in 1 cm steps)
Seat depth:	25 - 48 cm (in 1 cm steps)
Back height:	20 - 48 cm (in 1 cm steps)
Lower leg length:	20 - 50 cm (in 1 cm steps)
Rear seat height:	42, 45 cm
Front seat height:	44 - 54 cm (in 2 cm steps)
Tipping point:	6 - 10 cm (in 1 cm steps, option-dependent)
Total product weight:	Mechanical starting at 19 kg I electric starting at 35 kg
Max. standard load weight:	Up to 120 kg
Max. load weight for Heavy Duty:	Up to 140 kg appropriate individual dimensions upon request

FOLLOW US ON SOCIAL MEDIA









PRO ACTIV Reha-Technik GmbH

Im Hofstätt 11 72359 Dotternhausen info@proactiv-gmbh.de

www.proactiv-gmbh.com

Rev. 09/2024_01