



## YOUR ALL-ROUNDER FOR EVERY SITUATION IN LIFE.

The **SPEEDY 4all** active wheelchair offers you the perfect combination of sporty design and well-thought-out functionality. With its sturdy rigid frame, it stands for reliability and durability. Thanks to the integrated tipping point and seat height adjustment, you can flexibly adjust your sitting position at any time. Whether in everyday life or with special challenges – the **SPEEDY 4all** is ready to support you in every situation and master every challenge.

# SPEEDY 4all



### #MoreMobility

Traction devices and adaptive bikes can be connected in just a few steps via adapters; optional wheelbase extension for extra traction.



### #Reliability

Sturdy and durable frame for maximum mobility and compatibility with auxiliary drives and adaptations.



### #Independence

The flat frame design, the foldable backrest and the tool-free removable assemblies enable uncomplicated, independent transport.



### #Flexibility

Integrated tipping point and seat height adjustment in 1 cm steps for optimal adaptation to individual needs.







## PERFECTLY ADAPTED TO YOUR LIFE.

Your **SPEEDY 4all** is manufactured individually according to your needs. With 9 frame parameters, you can adapt your SPEEDY 4all perfectly to you. You have a large selection of footrests, brakes, push handles, back and seat systems, clothing guards, armrests, as well as drive and caster wheels.

### UNIQUE, SELECTABLE OPTIONS FROM PRO ACTIV:



The 7-way angle-adjustable and foldable backrest offers maximum flexibility for the ideal sitting position.



The footrest, which can be folded towards the back with a spring locking mechanism, is easy to handle and flexible in everyday life.



Manoeuvrable one-hand drive with double handrim for users with hemiplegia or restricted mobility on one side of the body.



With the Heavy Duty equipment, your SPEEDY 4all is designed for extreme loads and is suitable for user weights up to 190 kg.

## TECHNICAL SPECIFICATIONS

Seat width:	30 - 52 cm (in 1 cm steps)
Seat depth:	36 - 48 cm (in 1 cm steps)
Back height:	20 - 48 cm (in 1 cm steps)
Lower leg length:	37 - 50 cm (in 1 cm steps)
Rear seat height:	39 - 47 cm (in 1 cm steps)
Front seat height:	44 - 52 cm (in 1 cm steps)
Tipping point:	-3 - 10 cm (in 1 cm steps)
Total product weight:	Starting at 9.9 kg
Max. standard load weight:	Up to 120 kg
Max. load weight for Heavy Duty:	Up to 190 kg appropriate individual dimensions upon request

Data subject to design changes.

## FOLLOW US ON SOCIAL MEDIA



# PROACTIV®

PRO ACTIV Reha-Technik GmbH

Im Hofstätt 11  
72359 Dotternhausen  
info@proactiv-gmbh.de

[www.proactiv-gmbh.com](http://www.proactiv-gmbh.com)

Rev. 09/2024

Subject to technical modifications, printing errors & colour deviations.