



**Quotation date:** \_\_\_\_\_  
**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_  
**PRO ACTIV field representative:** \_\_\_\_\_

**Order date:** \_\_\_\_\_  
 according to GTC PRO ACTIV Reha-Technik GmbH

**Medical suppliers/company address:**

Cust. no.: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Street: \_\_\_\_\_  
 Postcode / city: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Fax: \_\_\_\_\_  
 E-Mail: \_\_\_\_\_  
 Contact: \_\_\_\_\_

**Shipping address:**  **same as company address**  
 **Other (see below)**

Name: \_\_\_\_\_  
 Street: \_\_\_\_\_  
 Postcode / city: \_\_\_\_\_

**Customer name/commission:**

Surname, forename: \_\_\_\_\_  
 Street: \_\_\_\_\_  
 Postcode / city: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Fax: \_\_\_\_\_  
 E-Mail: \_\_\_\_\_  
 D.O.B: \_\_\_\_\_  
 Paying organisation: \_\_\_\_\_  
 Indication: \_\_\_\_\_

**PRO ACTIV**  
**Transport restraint systems**

For the wheelchairs that have been tested and approved by PRO ACTIV

**Invoicing address:**  **same as company address**  
 **other (see below)**

Cust. no.: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Street: \_\_\_\_\_  
 Postcode / city: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Fax: \_\_\_\_\_  
 E-Mail: \_\_\_\_\_  
 Contact: \_\_\_\_\_


**Ordering the transport restraint systems for the following product:**

**New wheelchair order**  
 (See separate enclosed order sheet:  
 Commission: \_\_\_\_\_)


**Retrofitting to an existing PRO ACTIV wheelchair**

- Type: \_\_\_\_\_
- Product number: \_\_\_\_\_
- Year of constr.: \_\_\_\_\_
- Seat width SB: \_\_\_\_\_
- Seat depth ST: \_\_\_\_\_
- Back height RH: \_\_\_\_\_

## Dimension hip strap



Item no.	Description	Length A	Length B	Length C
E8000900007 A9000900001	Hip strap 1350 mm with lock	660 - 1410 mm	360 mm	230 mm
E8000900009 A9000900007	Hip strap 1670 mm with lock	910 - 1730 mm	480 mm	360 mm



Variant	Models affected	Dimension X
1	SPEEDY A1 / A2, TRAVELER, BUDDY classic / 4all	Approx. 30 mm
2	TRAVELER 4all Ergo, SPEEDY 4all / 4all Ergo, SPEEDY 4teen	Approx. 50 mm
3	SPEEDY 4you / 4you Ergo and TRAVELER 4you Ergo in the higher frame version	Approx. 90 mm

## Transport restraint systems for new wheelchair orders

### Transport restraint connection integrated on the rear frame with belt eye

€ 337.90

Frame-welded part for mounting a lap belt and the rear retractor. Lap belt and buckle (A9000900001 / A9000900007), please select separately if required.

- SPEEDY A1 / A2 TRAVELER A9000900030-K
- BUDDY classic A9040900001
- BUDDY 4all A9047900001



### Transport restraint connection on the rear frame via frame insert with belt eye For model: SPEEDY 4all / 4you / 4all Ergo / 4you Ergo, 4teen, TRAVELER 4all Ergo / 4you Ergo

€ 353.50

for mounting a lap belt and the rear retractor. Lap belt and buckle (A9000900001 / A9000900007), please select separately if required.

- A9000900020



## Accessories for new transport restraint orders (for all basic items - see above)


### Bearing block intermediate plate

Please note: The left wheel track guide increases by 12 mm per side.

with belt eye for mounting the front retractor

- Colour: black A9000900008-001 € 168.10
- Colour: Silver A9000900008-002 € 168.10



<p><b>1350 mm belt with lock, connection of transport restraints</b></p> <p><input type="checkbox"/> A9000900001      € 156.40</p> <p>Strap length adjustable. With belt links on both sides for mounting the diagonal shoulder belt. Belt colour: Red</p>	<p><b>1670 mm belt with lock, connection of transport restraints</b></p> <p><input type="checkbox"/> A9000900007      € 156.40</p> <p>Strap length adjustable. With belt links on both sides for mounting the diagonal shoulder belt. For use on wheelchairs with suitable seat shell. Belt colour: Red</p>	<p><b>Automatic shoulder belt with lock head</b></p> <p><input type="checkbox"/> A9000900002      € 129.10</p> <p>and deflector, length 3200 mm</p>
<p><b>Retractors for transport restraints</b></p> <p><input type="checkbox"/> A9000900003      € 496.10</p> <p>Connection to the wheelchair with belt loop or lock link. Connection to the vehicle floor via 4-way fitting.</p> <p>Set comprising: 2 retractors front: semi-automatic without hand tensioning wheel 2 retractors rear: semi-automatic with hand tensioning wheel</p>	<p><b>Headrest with segment padding</b></p> <p><input type="checkbox"/> A9000900004      € 384.30</p> <p>3-part and support tube ø 25 mm</p> 	<p><b>Headrest fixture</b></p> <p>for SPEEDY, BUDDY</p> <p><input type="checkbox"/> A9000900005-K001      € 258.60</p> <p>to accommodate the headrest with segment padding on the back cross bar. Reinforced version with double cross-tube.</p> <p>Detachable for TRAVELER</p> <p><input type="checkbox"/> A9000900005-K002      € 258.60</p> <p>to accommodate the headrest with segment padding on the end of the push handles. Reinforced version with double cross-tube.</p> <p><b>Please note:</b> for folding the wheelchair, the double cross-tube with headrest support must be unscrewed using the clamp lever.</p>

## Transport restraint systems for wheelchair retrofits

### Retrofit transport restraint connection

For model: **SPEEDY A1 / A2, TRAVELER, BUDDY classic**

Requirement for retrofitting: Wheelchair user less than 3 years  
Retrofitting is not possible if this requirement is not possible.

Integrated on the rear frame with belt eye.  
Frame-welded part for mounting a lap belt and the rear retractor.

R8000900004      € 954.60



For the retrofitting of the frame weld-part, it is necessary to carry out the following steps at our location:

- disassembly of the wheelchair
- decoating of the wheelchair frame
- welding on the frame welded part
- recoating the wheelchair frame
- assembly of the wheelchair

### Retrofit transport restraint connection



For model: **SPEEDY 4all / 4you / 4all Ergo / 4you Ergo / 4teen, TRAVELER 4all Ergo / 4you Ergo**

On the rear frame via frame insert with belt eye.

E8000900005      € 248.90

For mounting a hip strap as well as the rear retractor.  
Without assembly.



<b>Retrofit transport restraint connection</b> <b>For model: BUDDY 4all</b> with belt eye. <input type="checkbox"/> E8047900001 <b>€ 219.00</b> for connection to the rear retractor and for bolting to the lap belt. 		
<b>Accessories for transport restraint system retrofitting</b>		
<b>Bearing block intermediate plate</b> Please note: The left wheel track guide increases by 12 mm per side. with belt eye for mounting the front retractor and fixing screws <input type="checkbox"/> Colour: black E8000900046-001 <b>€ 266.00</b>		<input type="checkbox"/> Colour: Silver E8000900046-002 <b>€ 266.00</b> 
<b>1350 mm belt with lock, connection of transport restraints</b> <input type="checkbox"/> E8000900007 <b>€ 234.00</b> Strap length adjustable. With belt links on both sides for mounting the diagonal shoulder belt. Belt colour: Red	<b>1670 mm belt with lock, connection of transport restraints</b> <input type="checkbox"/> E8000900009 <b>€ 234.00</b> Strap length adjustable. With belt links on both sides for mounting the diagonal shoulder belt. For use on wheelchairs with suitable seat shell. Belt colour: Red	<b>Automatic shoulder belt with lock head</b> <input type="checkbox"/> E8000900003 <b>€ 201.90</b> and deflector, length 3200 mm
<b>Retractors for transport restraints</b> <input type="checkbox"/> E8000900014 <b>€ 772.30</b> Connection to the wheelchair with belt loop or lock link. Connection to the vehicle floor via 4-way fitting. Set comprising: 2 retractors front: semi-automatic without hand tensioning wheel 2 retractors rear: semi-automatic with hand tensioning wheel	<b>Headrest with segment padding</b> <input type="checkbox"/> E8000900004 <b>€ 598.10</b> 3-part	<b>Headrest fixture</b> for SPEEDY, BUDDY <input type="checkbox"/> E8000900006-K001 <b>€ 402.30</b> to accommodate the headrest with segment padding on the back cross bar. Reinforced version with double cross-tube. Detachable for TRAVELER <input type="checkbox"/> E8000900006-K002 <b>€ 402.30</b> to accommodate the headrest with segment padding on the end of the push handles. Reinforced version with double cross-tube. <b>Please note:</b> for folding the wheelchair, the double cross-tube with headrest support must be unscrewed using the clamp lever.
<b>Assembly of the transport restraint connection accessory by PRO ACTIV</b> <input type="checkbox"/> R8000900003 <b>€ 259.30</b>		

**Location, Date**

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

**Signature End customer/dealer**