

Im Hofstätt 11
D-72359 Dotternhausen
Phone +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
e-mail info@proactiv-gmbh.de
www.proactiv-gmbh.com



Quotation date: _____

**Consultation was carried out by
Staff at the medical supplies shop:** _____

PRO ACTIV field representative: _____

Order date: _____

According to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers/company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

Invoicing address, other(see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Shipping address, other (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name/commission:

Surname,
forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.O.B: _____

Paying
organisation: _____

Indication: _____



NJ1 adaptive bike

Up to 120 kg trailer load and 10 kg load capacity

- docking plate for quick coupling to the "Quick-Snap" adapter
- A-shaped handles with integrated switching and brake fittings
- secondary V-brake (rim brake) with locking mechanism
- parking stand
- chain protection tube
- black chain guard
- mudguard
- bottle holder with drinking bottle
- lighting set Sigma with rechargeable battery charger, front and rear lights

Please note: when purchasing an NJ1 adaptive bike in combination with a wheelchair, the reduced VAT rate applies for the entire unit.

Body size

Body height: _____ cm

Upper body
height: _____ cm

Weight: _____ kg

Information essential to process the order.

Drive unit

Black drive unit

- 20" (ETRTO 406) A9110100001-020 € 3,820.00
- 24" (ETRTO 507) A9110100001-024 € 3,820.00

Hub gear with integrated back-pedal brake

Shimano "Nexus" 8-speed hub gear

- A9110200023
307% transmission bandwidth
"Grip Shift" twist grip shifter
Gear ratio: Sprocket 32 teeth,
Cassette 20 teeth
(with bottom bracket support L – 19 teeth)



Dérailleur with disk brake

Shimano Deore 1x11 Rapid Fire 11-speed dérailleur

- A9110200032
mechanical dérailleur gear with operation on
the handle
11-speed cassette 11-42 teeth

Chain ring selection:

- 32 teeth recommended for
24" drive wheel
- 38 teeth recommended for
20" drive wheel
- 42 teeth for lower crank
frequency or driving
at higher speeds



Shimano Deore XT 1x11 Rapid Fire 11-speed dérailleur

- A9110200034 € 182.90
Mechanical dérailleur gear with operation on
the handle. 11-speed cassette 11-42 teeth
top chain ring with chain guards

Chain ring selection:

- 32 teeth recommended for
24" drive wheel
- 38 teeth recommended for
20" drive wheel
- 42 teeth for lower crank
frequency or driving
at higher speeds

Shimano Deore XT 1x11 DI2 with button operation 11-speed dérailleur

- A9110200035 € 1,414.50
electronic dérailleur with operation on the
handle. 11-speed cassette 11-42 teeth

Chain ring selection:

- 32 teeth recommended for
24" drive wheel
- 38 teeth recommended for
20" drive wheel
- 42 teeth for lower crank
frequency or driving
at higher speeds



Hub gear with disk brake

Shimano "Alfine" 11-speed hub gear

- A9100200026 € 563.50
"Rapid Fire" handle
Gear ratio: Sprocket 32 teeth,
Cassette 19 teeth



Rohloff "Speedhub 500/14" 14-speed hub gear

- A9100200027 € 1,426.00
"Grip Shift" twist grip shifter
Gear ratio: Sprocket 32 teeth,
Cassette 16 teeth



Gear change accessories

Mountain drive 2.5: 1 gear reduction for hills

A9100200003 € 759.00

Planetary gear set for doubling of the number of gears permits switching to the gear reduction for hills when stationary.
With 38 teeth chain wheel as standard.
Only in combination with grip size 40, 45, 48 and 53 cm.



Brake system

V-brake parking brake operated via locking lever

A9110300002 € 80.50

mounted on the bottom bracket support, brake handle on the handle, left side is omitted



Powerful brake system with brake disc Ø 203 mm

A9110300105 € 79.10

Recommended from a towed weight of 80 kg or body weight.
Not required with hydraulic back pedal disk brake selected
Not in conjunction with a Rohloff gear system



Double powerful brake system disc brake Ø 203 mm

A9110300100 € 330.40

including V-brake rim brake as parking brake, can be operated via locking lever
Recommended from a towed weight of 100 kg

Mechanical backpedal disk brake with crank release function

A9110300004 € 702.80

Compatible with dérailleur and free-running shifting hub.
In combination with:
- Grip width 45 cm, 50 cm, 53 cm or 58 cm
- Pedal crank crank lengths 155 or 175 mm
up to 120 kg body weight
Not in combination with Mountaindrive 2.5:1 gear reduction for hills

Tyres

20x1.75" Schwalbe "Marathon"

A9110105121
(ETRTO 406-47)
Touring tyres



20x2" Schwalbe "Knobby"

*(Low availability – only as long as stocks last.
Alternative option: Schwalbe "Marathon")*

A9110105126
(ETRTO 406-54)
Profiled tyres



Schwalbe "Marathon Plus"

20x1.75" (406-47) **€ 27.40**
A9130105128

24x1.75" (507-47) **€ 27.40**
A9110105135

unpuncturable, with anti-puncture belt
Smart Guard



Dimensions

Crank type and grip width

Bottom bracket shaft with Standard crank

Grip width 40 cm
A9110500100-040

Grip width 45 cm
A9110500100-045

Grip width 50 cm
A9110500100-050

Bottom bracket shaft with Tetra pedal crank (flared)

With flared crank arms for max. space for upper body

Grip width 48 cm (176 mm)
A9110500100-048

Grip width 53 cm (226 mm)
A9110500100-053

Grip width 58 cm (276 mm)
A9110500100-058

Standard pedal crank

Length 155 mm
A9110600100-155

Length 175 mm
A9110600100-175

Length 195 mm
A9110600100-195



Tetra pedal crank

Length 155 mm **€ 129.50**
A9100600105-155

Length 175 mm **€ 129.50**
A9100600105-175

Individual length 195 mm to 140 mm **€ 246.30**
A9100600100



Fixed bottom bracket support

L (X/Y 158/362 mm) **€ 63.50**
A9110501100-002

XL (X/Y 158/410 mm) **€ 63.50**
A9110501100-006

XXL (X/Y 170/465 mm) **€ 63.50**
A9110501100-007



Adjustable bottom bracket support

the angle and length can be adjusted via the clamped rotary point on the fork bridge and length.

For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.

M (300 mm) **€ 220.90**
A9110500001-001

L (370 mm) **€ 220.90**
A9110500001-002

XL (420 mm) **€ 220.90**
A9110500001-003

If no details are provided, the adjustable bottom bracket support shall be supplied in XL.

Deliverable up to size XL.



Handles

Palm rest on the handle

A9110601200 € 97.80

Padded, bolted, angle-adjustable.
Not in combination with other handles.



Tetra hand support, pivotable, adjustable grip distance

A9110601101 € 161.30



Tetra hand support, pivotable, adjustable grip width

A9110601102 € 182.40

Very stiff outer support as round tube
ø 20 mm with foam rubber cover
Cannot be pivoted in conjunction with Tetra
hand support



Tetra Ergo handle, free-standing

A9100600001 € 156.70

Crank handles arranged vertically
Can be combined with Tetra hand rest,
pivotable



Tetra shell handles

A9110601100 € 164.10



Bike accessories

Black aluminium mudguard

A9116500005 € 107.00

Only in combination with 20" drive wheel

Mudguard and chain guard unique colour selection

A9116500010 € 236.70

Made from a solid colour palette
Powder paint wheelchair

A9000100170-__ __ __

Colour: _____

Only in combination with 20" drive wheel

Bracket ø 22 mm x 42 mm for fixing of accessories (speedometer, lighting or similar)

A9110903271 Item price € 22.10

Colour: Natural anodised

_____ Items



Manoeuvring rollers for parking stands

A9110900017 € 116.20

Manoeuvring rollers mounted to the parking
stand via inlet



Manoeuvring bar mounted on the chain guard

A9110500000 € 160.80

with foam rubber cover



Aluminium luggage rack "cargo systems"

A9110903290 € 179.70

Colour: black
Max. load weight 10 kg on the
loading platform,
Max. 15 kg in combination with luggage bags
with low loading height.



<p>SIGMA BC 9.16 ATS Bike Computer</p> <p><input type="checkbox"/> A9110901000 € 66.80</p> <p>wireless, with bracket</p> 	<p>SIGMA BC 14.16 STS CAD Bike Computer</p> <p><input type="checkbox"/> A9110900018 € 163.20</p> <p>with pedal frequency set wireless</p> 
--	---

Options for users with reduced function

<p>Quick release clip lever extension for fixed frame adapter</p> <p><input type="checkbox"/> A9110700004 € 98.90</p> <p>to make the opening and closing of the quick release clips easier. Only in combination with rigid frame adapter.</p> 	<p>Securing the lever extension with locking bolt</p> <p><input type="checkbox"/> A9199000084 € 82.90</p> <p>Can only be used for the lever extension on the quick release clip.</p>	<p>Eccentric connection lever extension for Quick Snap connection simplifies opening or closing the Quick Snap Adapter. Tetra version</p> <p><input type="checkbox"/> A9110700005 € 70.30</p> <p>Straight version Only in combination with rigid frame adapter</p> <p><input type="checkbox"/> A9110700006 € 70.30</p> <p>Curved version Only in conjunction with front adapter</p> 
<p>Centring column adaptation assistance</p> <p><input type="checkbox"/> A9110700002 € 52.40</p> <p>in order to make inserting the central adapter easier.</p> 	<p>Parking stand with operating lever for high operating comfort</p> <p><input type="checkbox"/> A9110900004 € 55.90</p> <p>with longer actuating lever and reduced spring force.</p>	<p>Chest strap for upper body fixation</p> <p><input type="checkbox"/> Size S (105 cm) A9000800016 € 220.80</p> <p><input type="checkbox"/> Size M (125 cm) A9000800013 € 220.80</p> <p><input type="checkbox"/> Size L (145 cm) A9000800014 € 220.80</p> <p><input type="checkbox"/> Size XL (165 cm) A9000800017 € 220.80</p> <p>Consisting of elastic webbing material in combination with neoprene with Velcro and operating tab. Colour: black</p> 

Frame adapters

<p>SPEEDY fixed frame adapter with "Quick Snap Adaptation"</p> <p><input type="checkbox"/> A9110700141 € 492.50</p> <p>with quick release clips and centring column mount Colour: black</p> 	<p>TRAVELER folding frame adapter with "Quick Snap Adaptation"</p> <p><input type="checkbox"/> A9110701300 € 526.70</p> <p>with quick release clips and centring column mount on the folding mechanism Colour: black</p> 	<p>BUDDY fixed frame adapter with "Quick Snap Adaptation"</p> <p><input type="checkbox"/> A9110701400 € 492.50</p> <p>with quick release clips and centring column mount on the frame Colour: black</p> 
---	--	---

Wheelchair front adapter with universal clamping and "Quick Snap Adaptation"

A9110701500 € 561.10

Adapter mount on the wheelchair front end
Operation on the right side

Installation situation with PA products:
Frame stem G-shape from SB 31 cm to 54 cm
Frame stem V-shape from SB 30 cm to 50 cm

PLEASE NOTE:

If the tube diameter is not specified, as standard we shall deliver an adapter for a tube diameter of 23 - 27 mm.



Please make the following selection:

A9110701525

Tube diameter 23-27 mm

A9110701530

Tube diameter 28-32 mm

A9110701535

Tube diameter 33-37 mm

Clamp lever with 3-way star handle for actuating the frame clamping

A9199701507 € 101.00

mounted on the right and left for easy terminal operation with limited manual function
Only in conjunction with front adapter



Wheelchair-dependent accessories

Extended bike plate

For SPEEDY A1/A2 A9110100100-002 € 352.20

For BUDDY A9110100100-003 € 352.20

For TRAVELER A9110100100-001 € 352.20

4° camber integrated
Track can be adjusted independently



Wheelbase extension via "grid mechanism with pulling operation"

For SPEEDY 4all / 4teen / 4you / TRAVELER 4you A9110102001 € 537.50

removable, with continuous axle tube
4°/ 6° camber integrated
Track can be adjusted independently



Wheelbase extension, angle-adjustable and removable

For SPEEDY F2/F4 A9015100005 € 677.30

with continuous axle tube,
height adjustment via gear joint
adjustable in 6° angular increments
wheel position displacement approx. 14 cm
4°/6° camber integrated, adjustable track
Attachment possibility must be checked by PRO ACTIV.

Only in combination with base plate.
Up to 120 kg load capacity.



Services

Handover and instruction for the handbike/traction device by the dealer incl. mounting the frame adapter

on a PRO ACTIV wheelchair

A9990000008

Prod. no.: _____

on a third-party product

A9990000005

Manufacturer: _____

Model: _____

Frame tube ø ____ mm

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 7 - 2023-02-16 / The prices are plus statutory value-added tax and shipping costs = Included in the basic price = Position required
Order sheet NJ1 adaptive bike Rev. 05 UK