



Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____

Order date: _____
according to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: same as company address
 Other (see below)

Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:

Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.O.B: _____
Paying organisation: _____
Indication: _____

PRO ACTIV Transport restraint systems

For the wheelchairs that have been tested and approved by PRO ACTIV

Invoicing address: same as company address
 other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Ordering the transport restraint system for the following product:

New wheelchair order


(See separate order sheet enclosed:

commission.: _____)


Retrofitting to an existing PRO ACTIV wheelchair

- Type: _____
- Product number: _____
- Year of manufacture: _____
- Seat width SB: _____
- Seat depth ST: _____
- Back height RH: _____

Dimension lap belt



Article no.	Description	Lenght A	Lenght B	Lenght C
E8000900007 A9000900001	lap belt 1350 mm with lock	660 - 1410 mm	360 mm	230 mm
E8000900009 A9000900007	lap belt 1670 mm with lock	910 - 1730 mm	480 mm	360 mm



Variant	Models affected	Dimension X
1	SPEEDY A1 / A2, TRAVELER, BUDDY classic / 4all	approx. 30 mm
2	TRAVELER 4all Ergo, SPEEDY 4all / 4all Ergo, SPEEDY 4teen	approx. 50 mm
3	SPEEDY 4you / 4you Ergo and TRAVELER 4you Ergo in the version with high frame	approx. 90 mm

Transport restraint system for new wheelchair order

Transport restraint connection on rear frame integrated with belt eye

€ 354.80

Frame weld-on part for mounting the lap belt and for holding the rear retractor.
If required, please select lap belt with lock (A9000900001 / A9000900007) separately.

- SPEEDY A1 / A2 TRAVELER A9000900030-K
- BUDDY classic A9040900001
- BUDDY 4all A9047900001



Transport restraint connection on the rear frame via frame inlet with belt eye For models: SPEEDY 4all / 4you / 4all Ergo / 4you Ergo, 4teen, TRAVELER 4all Ergo / 4you Ergo

€ 371.20

for mounting the lap belt and for holding the rear retractor.
If required, please select lap belt with lock (A9000900001 / A9000900007) separately.

- A9000900020



Accessories for new transport restraint orders (for all basic articles - see above)




Bearing block intermediate plate




with belt eye for holding the front retractor.

Please note: The caster wheel track becomes 12 mm wider on each side.

- Colour: Black A9000900008-001 € 176.50
- Colour: Silver A9000900008-002 € 176.50



<p>Lap belt 1350 mm with lock, connection to transport restraint system</p> <p><input type="checkbox"/> A9000900001 € 164.20</p> <p>Adjustable belt length. With belt links on both sides for mounting the diagonal shoulder belt. Belt colour: Red</p>	<p>Lap belt 1670 mm with lock, connection to transport restraint system</p> <p><input type="checkbox"/> A9000900007 € 164.20</p> <p>Adjustable belt length. With belt links on both sides for mounting the diagonal shoulder belt. For use on wheelchairs with a seat shell. Belt colour: Red</p>	<p>Automatic diagonal shoulder belt with lock head</p> <p><input type="checkbox"/> A9000900002 € 135.60</p> <p>and deflection, length 3200 mm</p>
<p>Retractors for transport restraint system</p> <p><input type="checkbox"/> A9000900003 € 520.90</p> <p>Connection to the wheelchair with belt loop or lock link. Connection to the vehicle floor via 4-way fitting.</p> <p>Set consists of: 2 front retractors: semi-automatic without hand tensioning wheel 2 rear retractors: semi-automatic with hand tensioning wheel</p>	<p>Headrest with segment upholstery</p> <p><input type="checkbox"/> A9000900004 € 403.50</p> <p>3-part and support tube ø 25 mm</p> 	<p>Headrest holder ø 25 mm</p> <p>For SPEEDY, BUDDY</p> <p><input type="checkbox"/> A9000900005-K001 € 271.50</p> <p>to accommodate the headrest with segment padding on the back cross bar. Reinforced version with double cross-tube.</p> <p>Removable for TRAVELER</p> <p><input type="checkbox"/> A9000900005-K002 € 271.50</p> <p>to accommodate the headrest with segment padding on the end of the push handles. Reinforced version with double cross-tube.</p> <p>Please note: for folding the wheelchair, the double cross-tube with headrest support must be unscrewed using the clamp lever.</p> 
<p>Headrest holder with square mount for the support headrest</p> <p><input type="checkbox"/> Tube dimensions of support tube 15 x 15 mm € 281.40 A9000400850-001</p> <p><input type="checkbox"/> Tube dimensions of support tube 12,7 x 12,7 mm € 281.40 A9000400850-002</p> <p>Mounted on the back crossbar. Reinforced version with double cross tube. From SB 23 cm possible. Not in conjunction with central push handle.</p>		
<p>Transport restraint system for wheelchair retrofitting</p>		
<p>Retrofit for transport restraint system For models: SPEEDY A1 / A2, TRAVELER, BUDDY classic</p> <p>Requirement for retrofitting: wheelchair age under 3 years If this requirement is not met, a retrofit is not possible.</p> <p>Integrated on the rear frame with belt eye. Frame-welded part for mounting a lap belt and the rear retractor.</p> <p><input type="checkbox"/> R8000900004 € 1,002.30</p>  <p>For the retrofitting of the frame weld-part, it is necessary to carry out the following steps at our location:</p> <ul style="list-style-type: none"> - disassembly of the wheelchair - decoating of the wheelchair frame - welding on the frame welded part - recoating the wheelchair frame - assembly of the wheelchair 		

<p>Retrofit for transport restraint system For models: SPEEDY 4all / 4you / 4all Ergo / 4you Ergo / 4teen, TRAVELER 4all Ergo / 4you Ergo On the rear frame via frame inlet with belt eye.</p> <p><input type="checkbox"/> E8000900005 € 261.30 For holding the lap belt and the rear retractor. Without assembly.</p> 	<p>Retrofit for transport restraint system For model: BUDDY 4all with belt eye.</p> <p><input type="checkbox"/> E8047900001 € 372.00 for connection to the rear retractor and for mounting on the lap belt.</p> 	
Accessories for transport restraint system retrofitting		
<p>Bearing block intermediate plate with belt eye for holding the front retractor.</p> <p>Please note: The caster wheel track becomes 12 mm wider on each side.</p> <p><input type="checkbox"/> Colour: Black E8000900046-001 € 279.30</p> <p><input type="checkbox"/> Colour: Silver E8000900046-002 € 279.30</p>		
<p>Lap belt 1350 mm with lock, connection to transport restraint system</p> <p><input type="checkbox"/> E8000900007 € 245.70 Adjustable belt length. With belt links on both sides for mounting the diagonal shoulder belt. Belt colour: Red</p>	<p>Lap belt 1670 mm with lock, connection to transport restraint system</p> <p><input type="checkbox"/> E8000900009 € 245.70 Adjustable belt length. With belt links on both sides for mounting the diagonal shoulder belt. For use on wheelchairs with a seat shell Belt colour: Red</p>	<p>Automatic diagonal shoulder belt with lock head</p> <p><input type="checkbox"/> E8000900003 € 212.00 and deflection, length 3200 mm</p>
<p>Retractors for transport restraint system</p> <p><input type="checkbox"/> E8000900014 € 810.90 Connection to the wheelchair with belt loop or lock link. Connection to the vehicle floor via 4-way fitting.</p> <p>Set consists of: 2 front retractors: semi-automatic without hand tensioning wheel 2 rear retractors: semi-automatic with hand tensioning wheel</p>	<p>Headrest with segment upholstery</p> <p><input type="checkbox"/> E8000900004 € 628.00 3-part</p>	<p>Headrest holder ø 25 mm</p> <p>For SPEEDY, BUDDY</p> <p><input type="checkbox"/> E8000900006-K001 € 422.40 to accommodate the headrest with segment padding on the back cross bar. Reinforced version with double cross-tube.</p> <p>Removable for TRAVELER</p> <p><input type="checkbox"/> E8000900006-K002 € 422.40 to accommodate the headrest with segment padding on the end of the push handles. Reinforced version with double cross-tube.</p> <p>Please note: for folding the wheelchair, the double cross-tube with headrest support must be unscrewed using the clamp lever.</p>
<p>Assembly of the transport restraint connection accessory by PRO ACTIV</p> <p><input type="checkbox"/> R8000900003 € 272.30</p>		

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

Signature End customer/dealer