



FLEXIBLE, AGILE & READY FOR EVERY CHALLENGE.

SPEEDY 4teen

Distinct frame design, targeted weight optimisation and high-quality workmanship give the **SPEEDY 4teen** wheelchair for young people a sporty character. With 8 selectable frame parameters, it can be perfectly tailored to you. The SPEEDY 4teen easily adapts to your growing needs, thanks to the adjusting properties of the seat height and the tipping point. In this way, it supports you in your progress in dealing with the wheelchair and promotes your independence and mobility.




#GrowAlongFeature

The seat width increases from 26 to 38 cm by replacing assemblies, which makes the SPEEDY 4teen your reliable companion for many years.



#Reliability

The sturdy and durable frame for maximum mobility and compatibility with auxiliary drives and adaptations.



#MoreMobility

Traction devices and adaptive bikes devices can be connected in just a few steps via adapters.



#PerfectSittingPosition

You can choose between a straight or abducted frame. The abducted frame offers plenty of leg-room and allows for special leg positioning.



YOUR COMPANION FOR MORE MOBILITY.

The **SPEEDY 4teen** is specially manufactured according to your needs, both in terms of frame dimensions and features. You can choose from a variety of options, including footrests, brakes, push handles, back and seat systems, and drive and caster wheels. This makes your **SPEEDY 4teen** exactly the wheelchair that suits you perfectly.


UNIQUE, SELECTABLE OPTIONS FROM PRO ACTIV:



The 7-way angle-adjustable and foldable backrest offers maximum flexibility for the ideal sitting position.



The footrest, which can be folded towards the back with a spring locking mechanism, is easy to handle and flexible in everyday life.



Manoeuvrable one-hand drive with double handrim for users with hemiplegia or restricted mobility on one side of the body.



For handbike riding, more traction can be achieved via the wheelbase extension, which can be removed without tools.

TECHNICAL SPECIFICATIONS

Seat width:	26 - 43 cm (in 1 cm steps)
Seat depth:	32 - 48 cm (in 1 cm steps)
Back height:	20 - 48 cm (in 1 cm steps)
Lower leg length:	20 - 52 cm (in 1 cm steps)
Rear seat height:	36 - 47 cm (in 1 cm steps)
Front seat height:	40 - 50 cm (in 1 cm steps)
Tipping point:	1 - 10 cm (in 1 cm steps)
Total product weight:	Starting at 9.2 kg
Max. standard load weight:	Up to 100 kg

Data subject to design changes.

FOLLOW US ON SOCIAL MEDIA



PRO ACTIV[®]

PRO ACTIV Reha-Technik GmbH

Im Hofstätt 11
72359 Dotternhausen
info@proactiv-gmbh.de

www.proactiv-gmbh.com

Rev. 09/2024

Subject to technical modifications, printing errors & colour deviations.